

A vertical, clear plastic pipette tip is shown on the left side of the page. It has a grey cap at the top and a blue horizontal line indicating a liquid level. The tip is graduated with numbers 1, 2, 3, 4, and 5 from bottom to top. The background behind the tip is a light blue gradient.

# Important Information

Version B

## A new look for epT.I.P.S.<sup>®</sup> Macro Volume Tips

### **A new tip design**

Graduations on each tip allow easy visual recognition and control of liquid level when pipetting, dispensing or diluting with macro volume tips of 2.5 mL, 5 mL and 10 mL.

### **The new pipette tip**

The epT.I.P.S. 5 mL L with a length of 175 mm is made for easy reach of the bottom of e.g. 15 mL centrifugation tubes or graduated flasks, avoiding sample loss or cross contamination of pipette and vessel.

### **A new variant**

The epT.I.P.S. 5 mL pipette tips are now also available in Eppendorf purity grade Biopur<sup>®</sup>, free of human DNA and bacterial DNA, DNase-free, RNase-free, PCR inhibitor-free, ATP-free, Pyrogen-free and sterile for highest purity demands in cell culture, nucleic acid analytics, diagnostics or hygiene monitoring.

# New epT.I.P.S.<sup>®</sup> Macro Volume Tips Product Portfolio:

Volume range	Tip length	Purity grades	Packaging	ep Dualfilter T.I.P.S. <sup>®</sup> , PCR clean/Sterile
0.25–2.5 mL	115 mm	Eppendorf Quality™, PCR clean, Sterile (US only), Biopur	Standard (bulk), Reload, Box, Rack	No
0.1–5 mL	120 mm	Eppendorf Quality, Sterile (US only), Biopur <b>NEW</b>	Standard (bulk), Box, Rack	Yes
0.2–5 mL L <b>NEW</b>	175 mm	Eppendorf Quality, Biopur	Standard (bulk), Rack	Yes
0.5–10 mL	165 mm	Eppendorf Quality, Biopur	Standard (bulk), Rack	No
0.5–10 mL L	243 mm	Eppendorf Quality	Standard (bulk), Singles	Yes

**Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)**  
 Eppendorf AG · 22331 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com/consumables](http://www.eppendorf.com/consumables)